



Turbo Charge Your Microbial Lab

Introducing the 3M™ Petrifilm™ Plate Reader — Reads 3M™ Petrifilm™ Plates in Just 4 Seconds

3M Petrifilm Plates have already helped increase the productivity of QA personnel at more than 80 of the world's top 100 food companies by 82 percent. Now with the introduction of the new 3M Petrifilm Plate Reader — they can read and record results from those plates *in only 4 seconds*.

You get fast, automated readings of the three most common Petrifilm Plates — aerobic, coliform and e.coli/coliform counts. As a result, you not only increase productivity and reduce labor costs — but also eliminate variations in readings between lab techs.

Equally important, automated data storage saves time, eliminates transcription error, improves compliance and allows for easy export of data into Microsoft® Excel spreadsheets and uneditable log files.

For more information about how Petrifilm plates can help you turbo charge your lab, call us at 1-800-???-???. Or visit our website at www.3M.com/microbiology.



3M Microbiology